



# Selfist

REF D.19.50.59PK (div. Park's) - REF D.19.5059.B.PK (porta terza valva) **CE**

REF A.9023 - REF A.9045 - REF A.9060 - REF A.9080 **CE**  
0123

**Selfist** è una linea di valvette dentate, autostatiche, monouso e sterili, per fistulectomia. Sigillate singolarmente, sono disponibili in 4 lunghezze diverse (mm.23, 45, 60 e 80). Le **Selfist** vengono utilizzate in successione (dalla più corta alla più lunga) a mano a mano che la ferita diventa più profonda e garantiscono la visibilità ottimale del fondo ove procede la fistulectomia. In tal modo l'uso di divaricatori da parte di uno o due assistenti risulta superfluo.

Le **Selfist** possono essere innestate all'impugnatura in acciaio del **divaricatore Park's** e al **porta terza valva** marchiati Sapi Med o, in alternativa, ai divaricatori Park's pluriuso muniti di **attacco a pressione** delle valve.

**Latex free.**

**Conf.vendita REF D.19.50.59PK (divaricatore Park's): 1 pz.**

**Conf.vendita REF D.19.5059.B.PK (porta terza valva): 1 pz**

**Conf.vendita REF A.9023 (L. mm.23): scatola 24 pz.**

**Conf.vendita REF A.9045 (L. mm.45): scatola 24 pz.**

**Conf.vendita REF A.9060 (L. mm.60): scatola 24 pz.**

**Conf.vendita REF A.9080 (L. mm.80): scatola 24 pz.**





# Selfist

REF D.19.50.59PK (Park's retractor) - REF D.19.5059.B.PK (third blade holder)



REF A.9023 – REF A.9045 – REF A.9060 – REF A.9080



**Selfist** is a line of self-retaining, single-use, sterile, serrated blades, for fistulectomy. Sealed singly, available in 4 different sizes (23, 45, 60 and 80 mm). **Selfist** instruments are used in succession (from shortest to longest) as the wound gets deeper and guarantee an excellent view of the bottom where the fistulectomy is taking place. In this way, there is no need for one or more assistants to use retractors.

**Selfist** instruments can be fitted into the steel handle of the Park's retractor and into the third blade holder branded Sapi Med or, alternatively, onto Park's multi-use retractors equipped with pressure blade attachments.

**Latex free.**

**Pack REF D.19.50.59PK (Park's retractor): 1 piece.**

**Pack REF D.19.5059.B.PK (third blade holder): 1 piece**

**Pack REF A.9023 (L. 23 mm): box 24 pcs.**

**Pack REF A.9045 (L. 45 mm): box 24 pcs.**

**Pack REF A.9060 (L. 60 mm): box 24 pcs.**

**Pack REF A.9080 (L. 80 mm): box 24 pcs.**

